STATIS[®]

AM 20

Operating instructions
Manuel de l'utilisateur
Benutzerhandbuch
Manuale per operatori
Manual del Operador

Caution: In the U.S.A., Federal Law restricts this device to sale by or on the order of a dentist.

SciCan Dental

Your Infection Control Specialist™

Distributed by:

SciCan Ltd. 1440 Don Mills Road, Toronto ON M3B 3P9 CANADA

www.scican.com

Phone: +1 416 445 1600 Toll free: 1-800-667-7733 Fax: +1 416 445 2727

Manufactured by:

SciCan GmbH

Wangener Strasse 78

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GERMANY

Tel.: +49 (0)7561 98343 - 0

Fax: +49 (0)7561 98343 - 699

SciCan Medtech

Alpenstrasse 16

6300 Zug

SWITZERLAND

Phone: +41 (0) 41 727 7027

Fax: +41 (0) 41 727 7029

SciCan Inc.

701 Technology Drive

Canonsburg, PA 15317

USA

Phone: +1 724 820 1600

Fax: +1 724 820 1479

Toll free: 1-888-632-6031

Introduction

Dental air driven motor STAT/S™ AM 20 with internal spray, without light as per ISO Standard 3964. Use in dentistry for prophylaxis, general dentistry and endodontics work.

Important Information

Do not use the motor at pressures in excess of the maximum. The air pressure required for use is between 3.1 – 3.9 bars measured during operation of the device at the rear of the instrument.

Application

Product intended for professional use only.

The technical specifications, illustrations and dimensions contained in these instructions are given only as a guide. They may not be the subject of any claim. The manufacturer reserves the right to make technical improvements to its equipment, without amending these instructions. For all additional information, please contact SciCan at the address indicated on the front page.

The device must be used by a competent person in accordance with the current legal provisions concerning industrial safety, health and accident prevention regulations, and these working instructions. In accordance with these requirements, the operator:

- must only use working devices that are in perfect working order; in the event of irregular functioning, excessive vibration, abnormal heating or other signs indicating malfunction of the device, the work must be stopped immediately; in this case, contact a repair centre that is approved by SciCan;
- must ensure that the device is used only for the purpose for which it is intended, must protect himself, his patients and third parties from any danger, and must avoid contamination through the use of the product.

The device is intended for medical treatment only; any use other than that for which this product is intended is not allowed and may prove to be dangerous. The medical device meets all the current national legal requirements.

It is essential to use dry, purified compressed air in order to ensure the long working life of the device. Maintain the quality of the air and the water by regular maintenance of the compressor and the filtration systems. The use of unfiltered hard water will lead to early blockage of the tubes, connectors and spray cones.

Rest the device on a suitable support in order to eliminate the risk of injury.

Guarantee

The manufacturer grants the operator a warranty covering all functional defects, material or production faults. The device is covered by this warranty for 2 years from the date of invoicing. In the event of justified claims, the manufacturer or its authorised representative will fulfill the company's obligations under this guarantee by repairing or replacing the product free of charge. Any other claims, of whatever nature, in particular in the form of a claim for damages and interest, are excluded.

The manufacturer shall not be held responsible for damage or injury and the consequences thereof, resulting from:

- excessive wear and tear
- improper handling
- non-observance of the instructions for installation, operation and maintenance
- unusual chemical, electrical or electrolytic influences
- poor connections, whether of the air water or electricity supply.

The guarantee shall become null and void if the damage and its consequences are due to improper manipulation of the product, or modifications to the product carried out by persons not authorized by the manufacturer.

Claims under the terms of the guarantee will be considered only on presentation, together with the product, of the invoice or the consignment note, on which the date of purchase, the product reference and the Serial No. should be clearly indicated.

Installation

Place the bottom of the AM 20 motor into female end of dental units tubing and secure in place by rotating tubing collar onto the motor.

Maintenance

Non-sterilizable. When in use, use a protective sheath on the AM20 motor.

Only use maintenance products and/or components from SciCan i.e. SciCan STATCARE™ spray or the SciCan automatic lubrication device STATMATIC™. The use of other products and components can lead to annulment of the guarantee.

After cleaning and disinfecting the air motor and before using it, operate at moderate speed for 10 to 15 seconds to distribute, and remove, the excess lubricant.

Servicing

Never disassemble the device. For all modifications and repairs, we recommend that you contact your regular supplier or SciCan directly. The manufacturer recommends inspection of the instruments on a regular basis.

Recommended reprocessing methods according to DIN EN ISO 17664

The operator should be aware of personal health risks when reprocessing contaminated dental handpieces and air motors should wear gloves and eye protection.

- Wipe off any debris or blood contamination at chair side.
- Do not bathe dental handpieces and air motors in disinfection solutions.
- Reprocess immediately after use.
- Remove the attachment before reprocessing!
- NON-sterilizable!

Cleaning and Disinfection

Please note that the AM20 should not be cleaned in an ultrasonic bath.

Attention: Do not use products containing acetone, chlorine or bleach as disinfectants

Manual External Cleaning and Disinfection

Cleaning:

After each use, the AM20 must be thoroughly scrubbed either under water (potable water) or with 60 – 70% alcohol solution.

Disinfection

Not use products containing acetone, chlorine or bleach as disinfectants! When applying disinfection sprays/wipes, follow the instructions for use of these materials.

FOR GERMAN USERS:

In order to efficiently reprocess the AM20, internal cleaning and disinfection of the spray channels can be achieved through the use of ALPRO WL-clean*, ALPRO WL-dry* and ALPRO WL-cid*.

*ALPRO WL-clean, ALPRO WL-cid and ALPRO WL-dry are registered trademarks of ALPRO MEDICAL GMBH, Mooswiesenstrasse 9, 78112 St. Georgen, Germany.)

Drying

Manual Drying

Apply dry air to the internal and external surfaces until surface water is no longer visible.

Lubrication

After cleaning and / or disinfecting the the AM20 it is important to maintain/lubricate with SciCan Statcare Spray or use an automated care system such as the Statmatic maintenance unit.

Manual Lubrication

- Place the motor in a cloth or paper towel so as to catch any debris and foreign matter expelled by the atomizer jet.
- Insert the nozzle of the lubricant into the back of the instrument sleeve.
- Spray the STATCARE Spray for 1-2 seconds.
- Run the handpiece for 15-20 seconds to remove any residual debris and oil.

Sterilisation of the STATIS AM20

Please note that the STATIS AM20 motor CANNOT be sterilized!

Storage

Reprocessed products should be stored wrapped in a dust-free, dry, dark environment at room temperature.

Environment

Temperature between -40°C (-40°F) and 70°C (158°F), relative humidity between 10% and 100%, atmospheric pressure 50 kPa to 106 kPa (7,3 to 15,3 psi).

Specifications

Tubing:	4-hole tubing
Dimensions:	Ø 22.8mm Length 103.7
Weight:	84 gr
Mechanical power:	16 W
Torque:	2.4 Ncm at 0 rpm
Noise level:	less than 70 dBA at 45 cm in accordance with ISO 7785-1
Air consumption noutput:	56 NI

In these instructions, motor corresponds to the product described in the heading STAT/SAM 20.

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STAT/S, STATCARE, and STATMATIC are trademarks of SciCan Ltd. OPTIM is a registered trademark of SciCan Ltd.
STAT/S, STATCARE et STATMATIC sont des marques de fabrique de SciCan Ltd. OPTIM est une marque déposée de SciCan Ltd.
STAT/S, STATCARE und STATMATIC sind Marken von SciCan Ltd. OPTIM ist eine eingetragene Marke von SciCan Ltd.