



Think of it as a routine checkup for your sterilizer.

The ENSURE line of infection control monitoring products ensures your equipment is performing at its peak. It's a simple, easy to use process to validate your equipment. The unique incubator allows you to monitor up to 5 sterilizers at one time, providing you with easy to read results that show whether sterility has been achieved... or not. In addition, Class 4 chemical indicators are also available for your steam autoclave. The Ensure system can satisfy even the most stringent international regulatory requirements. No need to wonder, when it's this easy to know.

Time for a checkup? Visit www.scican.com

Your Infection Control Specialist™

SciCan Dental



EnsureTM ✓

Because the Stakes are too highTM

The Ensure product line is an easy to use, in-office monitoring system for steam sterilization.

CHEMICAL INDICATORS (CI)

The ENSURE CI is a fast, economical solution for load monitoring. SciCan offers a Class 4 indicator to meet the most stringent requirements in today's changing regulatory environment.



ENSURE CHEMICAL INDICATORS (CI)

Part #	Description
SCI134 Class 4	Chemical indicator strips 132-134°C (250 count)

BIOLOGICAL INDICATORS (BI)

The ENSURE BI's are fast, reliable, self-contained indicators for results in as little as 24 hours. BI's are the only cost effective method to ensure effective performance of your sterilizer.

The ENSURE Incubator, a 10-well incubator, is a very unique unit unlike any in the market today. With an integrated vial crusher and built in timer for up to five different programs, the SciCan ENSURE BI system leaves no room for error in time.

The ENSURE BI's are available in a starter kit system that consist of 25 BI's, the unique incubator and an easy to fill-in log book. The starter kit is the most cost effective, accurate system for monitoring your steam autoclave directly in your office in only 24 hours.



ENSURE BIOLOGICAL INDICATORS (BI)

Part #	Description
SCI-BI25	Ensure BI's (25 count refill box)
SCI-BI100	Ensure BI's (100 count refill box)
SCI-BISK-120	Ensure system, Starter Kit (25 BI's, incubator, log book)
SCI-BINC-120	Ensure incubator (incubator only)
SCI-BILOG	Ensure biological monitoring system log book

SciCan Inc., 701 Technology Drive, Canonsburg, PA 15317
 Phone (724) 820-1600; Fax (724) 820-1479; Toll Free 1-800-572-1211
www.scican.com