




STATFlow™

Never empty your condenser bottle again.

STATFlow, the only direct-to-drain system designed from the ground up for your STATIM® 2000 and 5000 cassette autoclaves.



STATFlow
Code: 01-113025

FAST
Drains STATIM exhaust safely and automatically.

EFFICIENT
Eliminates wasted staff time emptying condenser bottle.

COST EFFECTIVE
Dedicate more resources to patient care and less on maintenance.

FEATURES

UNIVERSAL SYSTEM. STATFlow will work with either the STATIM 2000 or 5000.

SMART DESIGN. Triangular water tank is shaped to take up less space under your sink.

ECO DESIGN. Water tank is made of recycled plastic material.

EASY TO USE. No consumable parts and no daily maintenance is required.

ROBUST. System will continuously cool your STATIM's hot exhaust and drain it safely – no matter how many cycles you run.

COMPLIANT. All components and system set-up are compliant with most regional plumbing codes.

VALIDATED. STATFlow is the only direct-to-drain product designed to safely and effectively drain the steam exhaust from a STATIM 2000 or 5000.

COMPATIBLE. Can be used with all drainage pipe material types.

Please consult with your authorized SciCan dealer or representative for pre-installation checklist.