



STATIM® SILICONE MAT

Keep your STATIM looking new with the STATIM Silicone Mat!



STATIM 2000 G4
Code: 01-112867S

The NEW STATIM Silicone Mat

provides a much safer and hygienic working area by securely seating the cassette on top of the STATIM. The mat also protects paintwork from accidental scuffs and scrapes, keeping your Statim looking new for longer. The high grade silicone also helps dissipate heat and protects surfaces from damage caused by heat transfer from hot metal surfaces.

Available for:
STATIM 2000, 5000 and 7000.



STATIM 5000 G4
Code: 01-112868S



Contact SciCan for more details.
www.scicanusa.com